


Evaluation CHART

Use this chart to determine your level of hair loss.

Shapiro Medical will help you determine the number of grafts needed per session.



The number of follicular units transplanted per session and the post-op density it yields

Determined by multiplying the surface area (number of sq cm) by number of follicular units we can transplant per sq cm, per session (20-25 fu/sq cm)

FU = Follicular Units

* Thinning areas 1 & 2 may require additional grafts for fullness
Note: results based on working on an area of 100 sq/cm (frontal area).

DEGREE OF THINNING SCALE					
	1*	2*	3	4	5
PRE-SURGERY DENSITY No. of FU/sq. cm.	 BALD 00 FU	 10-20 FU	 20-25 FU	 40-50 FU	 NO THINNING 80-100 FU
NUMBER OF GRAFTS NEEDED PER SESSION					
FRONTAL AREA					
Hairline zone (50 sq cm)	1,000-1,250	1,250	1,000	—	—
Frontal core (50 sq cm)	1,000-1,250	1000	1,000	—	—
TOTAL FRONT HALF (100 sq cm)	2,000-2,500	2,000-2,250	2,000		
POST SURGERY DENSITY No. of FU/sq. cm.	20-25 FU	30-45 FU	40-50 FU	—	—
CROWN					
Softball	1,500	1,200	—	—	—
Baseball	1,200	1,000	—	—	—
Golfball	800	500-800	—	—	—